



COOLTOUCH™ CTEV 1320™ Endovenous Laser

Nd:YAG Laser for Endovenous Closure of the Great & Small Saphenous Veins



SPECIFICATIONS

- Wavelength:** 1320 nm
- Power:** Up to 10 Watts
- Delivery:** 600µm, 273µm and 2.5F Fibers
- Dimensions:** 31" high x 18" wide x 21" deep
(79 cm x 46 cm x 53 cm)
- Weight:** 160 lbs (73 kg)
- Utilities:** 115 VAC, 50/60 Hz, Single Phase
230 VAC, 50/60 Hz, Single Phase

SELECT FEATURES

- 1320 nm wavelength targets water in vessel wall
- More comfortable treatment for patients
- Little to no post-op bruising, swelling or pain
- Unique automated fiber Pull-Back Device
- Aiming beam allows identification of treatment area

CoolTouch Pull-Back Device



9085 Foothills Blvd., Roseville, CA 95747
877.858.2665 (toll free) • 916.677.1900 (phone) • 916.677.1901 (fax) • www.cooltouch.com

